

# AQUAGLO FACIAL

## 3 STEP EXPRESS FACIAL

Cleanse, Exfoliate, Extract and Hydrate in an instant



### INDICATIONS:

All skin types

### SUGGESTED FREQUENCY OF TREATMENT:

1 x week for 6 treatments. 1 x month for maintenance.

### CLIENT PREPARATION:

Remove all jewellery and protect client's hair with a cap or headband.

Wash your hands and put gloves on.

### CLEANSING PROCEDURE:

Remove eye and lip make-up with an appropriate product  
Cleanse the face thoroughly with a foaming or gel style cleanser. Remove cleanser with sponges or a warm, aromatic towel. Apply Toner if required.

1

### STEP 1

#### DEEP CLEANSE & EXFOLIATE - HYDRADERMABRASION

Use the small spiral tip to exfoliate the skin with Solution 1 (Cleanse/Exfoliate solution).

Press and hold the tip on the back of your hand to bring solution to the surface.

Perform a thorough pass covering all areas of the face with a 50% overlap on each stroke. Avoid the orbital area. Use sweeping strokes similar to traditional Microdermabrasion.

Work from the centre areas to the outer areas to encourage Lymphatic Drainage.

Use very low pressure to include the neck area

Repeat process a second time on any problem areas.

*Keep small tip for Extraction (Step 2).*

2

### STEP 2

#### GENTLE ACID PEEL & EXTRACTION

Use the small tip and Solution 1 (Cleanse/Exfoliate Solution).

Concentrate on any areas where congestion or comedones are present. Use a higher vacuum setting, firmer stroking pressure and a slight "twist and press" movement of your hand to extract impurities from the skin.

*Remove small spiral tip and set aside for cleaning and*

*Sterilization, or dispose of it*

3

### STEP 3

#### CALMING & HYDRATING LIQUID INFUSION

Use large spiral tip to infuse Solution 3 (Calming).

Press and hold the tip on the back of your hand to bring solution to the surface.

Use a lower vacuum setting with a "wet kiss" (dab and press) motion to saturate skin with solution. Cover the entire face and neck with this movement (avoid orbital area). Use your fingers to lightly rub in any excess solution on the skin surface.

*Remove large spiral tip and set aside for cleaning and Sterilization, or dispose of it*

### FINAL HYDRATION

Apply small amount of Treatment Cream, Eye Cream and Sun Protection to finish.

# DELUXE GLOW FACIAL

INCORPORATED WITH YOUR TARGETED FACIAL PRODUCTS

*Cleanse, Exfoliate, Hydrate and Rejuvenate with a Deluxe Intensive Facial*



1

**STEP 1**  
**DEEP CLEANSE & EXFOLIATE**  
**- HYDRADERMABRASION**

Use the small spiral tip to exfoliate the skin with Solution 1 (Cleanse/Exfoliate solution).

Press and hold the tip on the back of your hand to bring solution to the surface.

Perform a thorough pass covering all areas of the face with a 50% overlap on each stroke. Avoid the orbital area. Use sweeping strokes similar to traditional Microdermabrasion. Work from the centre areas to the outer areas to encourage Lymphatic Drainage.

Use very low pressure to include the neck area  
Repeat process a second time on any problem areas.

*Keep small tip for Extraction (Step 2).*

2

**STEP 2**  
**GENTLE ACID PEEL & EXTRACTION**

Use the small tip and Solution 1 (Cleanse/Exfoliate Solution).

Concentrate on any areas where congestion or comedones are present. Use a higher vacuum setting, firmer stroking pressure and a slight "twist and press" movement of your hand to extract impurities from the skin.

*Remove small spiral tip and set aside for cleaning and Sterilization, or dispose of it*

3

**STEP 3**  
**CALMING & HYDRATING LIQUID INFUSION**

Use large spiral tip to infuse Solution 3 (Calming).

Press and hold the tip on the back of your hand to bring solution to the surface.

Use a lower vacuum setting with a "wet kiss" (dab and press) motion to saturate skin with solution. Cover the entire face and neck with this movement (avoid orbital area). Use your fingers to lightly rub in any excess solution on the skin surface.

*Remove large spiral tip and set aside for cleaning and Sterilization, or dispose of it*

4

**STEP 4**  
**MASSAGE**

Apply a small amount of Massage Product and perform Traditional massage, Lymphatic Drainage or Pressure Points on the face, neck and décolleté. Remove excess product with a warm, aromatic towel (optional).

5

**STEP 5**  
**MASK**

Apply a mask of your choice to infuse the Calming Solution S3 even further, and to intensively hydrate and restore the skin. Apply to the face and neck. Any type of Cream, Gel, Clay, Collagen Fleece, Rubber or setting masks may be used, depending on skin type and products available. Eye masks can be included as an extra treatment option.

*Remove mask after 15 minutes, using sponges or a warm aromatic towel.*

**FINAL HYDRATION**

Apply small amount of Treatment Cream, Eye Cream and Sun Protection to finish.